

Guide For FHS/T00 and FHS/T01



FHS/T00



FHS/T01



MG - 18; MG - 23; MG - 27



MAGNIFIER LENSES



MG - 18; MG - 23; MG - 27

EMISSION FILTER - FS/TEF-xxx

Emission Filter Model Number	Fluorescent Proteins	Open (Wavelength nm)
FS/TEF-1C2	Cyan	480-505 nm
FS/TEF-2G2	Green	500-515 nm
FS/TEF-3GY1	Green & Yellow	500-650 nm
FS/TEF-3GY2	Green & Yellow	500-550 nm
FS/TEF-4Y1	Yellow	520-650 nm
FS/TEF-4Y2	Yellow	520-550 nm
FS/TEF-4Y3	Yellow	520-540 nm
FS/TEF-4R2	Red	590-660 nm
FS/TEF-4O2	OrangeRed	567-590 nm



Camera/Video Stand
CVS 8

FHS/T00 can be used with:



Excitation stand
FBL/Basic-B&N-01



Flashlight
UFP-MDS-CC/HB

LIGHT SOURCE
FS/ULS-02-cc



Light Source Model Number	Excitation Wavelength	Emission for Fluorescent Proteins
FS/ULS-02RB	440-460 nm	Cyan
FS/ULS-02B2	460-495 nm	Green & Yellow
FS/ULS-02G2	515-555 nm	Red
FS/ULS-02OR	567-588 nm	OrangeRed



BIOLOGICAL LABORATORY EQUIPMENT
MAINTENANCE AND SERVICE LTD.

www.bls-ltd.com

Zsélyi Aladár u.31, H-1165 Budapest, HUNGARY

Tel: +36-1-407-2602 Fax: +36-1-401-0925 E-mail: bls@t-online.hu Website: www.bls-ltd.com